

Portable Suction Pump

EAR WAX REMOVAL BY MICRO-SUCTION









Light Weight & Compact

Suction Cannister

Bacterial Filter

Silicone Tubing

Portable Suction Pump

EAR WAX REMOVAL BY MICRO-SUCTION

The VMS-PK Suction Pump is a fully portable electrically powered device that is ideal for effective ear wax removal by micro-suction. It offers a gentle suction flow rate of 15-20L/min and a vacuum level of 500-600mmHg. The VMS-PK Suction Pump is light weight, quiet and fully portable.

Light Weight & Fully Portable

The VMS-PK Suction Pump is fully portable, weighs just 3.9kg and features compact dimensions of 350mm x 260mm x 290mm.

Quiet Operation

The VMS-PK operates quietly at less than 53 dB sound level.

Optimum Flow Rate

The VMS-PK Suction Pump operates between 15 – 20 L/min.

Gentle & Effective Vacuum Level

The VMS-PK Suction Pump can be set to operate between 0.065 – 0.080 Mpa (approx 500-600 mmHg) for effective micro-suction performance.

Simple Operation

The VMS-PK is simple to set up and maintain. Always read the Instructions For Use supplied.

Sturdy Design

The VMS-PK features anti over-flow technology; water proof switch; shatter - proof & reusable polycarbonate cannister.

Micro-suction Accessories

Vorotek stocks a wide range of accessories required for micro-suction including the VorOtek O Scope for direct vision, suction tips, controllers, tubing & spare parts.



Vorotek Pty Ltd 23 Point Henry Road MOOLAP VIC 3224 Australia

T +61 3 5248 5809 E office@vorotek.com.au www.vorotek.com

