



# ENDO Light

PORTABLE LED ILLUMINATION FOR ENDOSCOPES

Convenience without compromising Illumination

Made in Australia

**VOROTEK**  
experience. knowledge. performance



# ENDOLight

## Convenience without compromising Illumination

The Vorotek ENDOLight is replacing fibre optic lighting systems for endoscopes. The intensity & quality of the illumination is exceptional and not being tethered by a fibre optic cable or power lead is a major benefit.

### Surgical Quality LED Illumination

Vorotek's LED (50,000 hours operating life) produces a brilliant **150,000 lux**. The intensity and consistency of the illumination is comparable to fibre optic systems.

### Dimmable LED Illumination

The dimmer function provides 8 evenly spaced levels of illumination.

### Endoscope Compatibility

For flexible and rigid Endoscopes. Fits directly to **Pentax** endoscopes. Adaptors available for **Olympus** and **Storz** endoscopes.

### Interchangeable

The ENDOLight can be used with any Vorotek Battery Pack e.g. the *Vorotek O Scope Battery Pack*.

### Ease of Use

Portable, compact and light weight, emphasizing simplicity and durability.

### Portable Power Supply

The high-capacity Lithium Ion Power Pack (either belt or pocket worn), provides 7hrs to >20hrs "on time" depending on the illumination level.

### Durable

Outstanding reliability and durability results from experience and knowledge acquired by Vorotek over 40 years.

### Designed & Manufactured in Australia

Vorotek Pty Ltd  
23 Point Henry Road  
MOOLAP VIC 3224, Australia

T +61 3 5248 5809  
E [office@vorotek.com.au](mailto:office@vorotek.com.au)  
[www.vorotek.com](http://www.vorotek.com)

**VOROTEK**  
experience. knowledge. performance