

VOROTEK

MICRO-SUCTION & CURETTAGE START-UP KIT

Setting up a clinic with micro-suction equipment is simple and cost-effective. The VorOtek Micro-Suction & Curettage Start-Up Kit is perfect to support you in making the transition to micro-suction under direct vision or for setting up in a new clinic.

The start-up kit includes the following equipment:

VorOtek O Scope.

HeadBand or SpecFrame model. 12mm Converged Binocular Optics & Fully Integrated LED Illumination with 3 dioptre Magnification. Adjustable IPD. Includes Power Pack, Charger, Carry Case



Speculum Set.

15 speculums (3 of each size XS, S, M, L, XL). Autoclavable plastic



Wax Rings & Hook.

2.5mm & 3.5mm size wax rings
Small 90 degree hook



Suction Pump.

Flow 15-20L/min; max vacuum 600mmHg). Includes x1 suction tubing, silicone extension tubing & connector, micro ear suction tip (green), controller handle



Micro Ear Suction Tips.

Green 14 Gauge (2.1mm OD x 70mm) pack of 50 sterile, single use



Controller Handles.

Pack of 10 NON-Sterile



Start-Up Kit price: \$2,782* ex GST (order code VMS-SUP1)

Please contact VorOtek for all enquiries.

If you would like a quote or to place an order, please go to <https://www.vorotek.com/order-quote> select micro-suction in step 1 & complete the Order/Quote Form steps.

*All prices exclude GST, postage & handling & apply to Australian customers only.

Vorotek Pty Ltd
T: +61 3 5248 5809
E: office@vorotek.com.au
www.vorotek.com